


Product Information	
Catalog No.	BS11957
Product name	Anti-Human CD357/TNFRSF18/GITR Reference Antibody (EU102)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	EU102
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TL6, TNFSF18, Activation-inducible TNF-related ligand, GITRL, Glucocorticoid-induced TNF-related ligand, Tumor necrosis factor ligand superfamily member 18, hGITRL, AITRL
Expression system	Mammalian Cells
Accession	Q9Y5U5
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	EU102
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human CD357/TNFRSF18/GITR (EU102)

SDS-PAGE