


Product Information	
Catalog No.	BS12023
Product name	Anti-Human PLA2G1B/PLA2 Reference Antibody (DIACC 1010)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	DIACC 1010
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Phospholipase A2, PLA2G1B, Group IB phospholipase A2, PLA2A, PPLA2, Phosphatidylcholine 2-acylhydrolase 1B, PLA2
Expression system	Mammalian Cells
Accession	P04054
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	DIACC-1010, DIACC1010, DIACC1010
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human PLA2G1B/PLA2 (DIACC 1010)

SDS-PAGE