


Product Information	
Catalog No.	BS12039
Product name	Anti-Human IL12B/IL-12 p40/NKSF2 Reference Antibody (ART123)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	ART123
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL-12 subunit p40, IL12B, Interleukin-12 subunit beta, IL-12B, CLMF p40, NKSF2, NK cell stimulatory factor chain 2, Cytotoxic lymphocyte maturation factor 40 kDa subunit
Expression system	Mammalian Cells
Accession	P29460
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	ART123
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human IL12B/IL-12 p40/NKSF2 (ART123)

SDS-PAGE