

Product Information

Catalog No.	BS12132
Product name	Research Grade Anpocogin

Background

Host species:

Isotype: Ancylostomacaniumnematodeanticoagulantproteinc2.

Species reactivity: Ancylostomacanium(Doghookworm)

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality: Monoclonal

Applications: Research Grade Biosimilar

Target: AnticoagulantproteinC2

Purification: Purified by Nickel column.

Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: CAS:2725767-44-0

Note: For research use only. Not suitable for clinical or therapeutic use.