

Product Information

Catalog No.	BS12133
Product name	Research Grade Berahyaluronidase Alfa

Background

Host species:

Isotype: Human hyaluronidase isozyme-5/HYAL5 fragment engineered variant with the structural domain replaced by the corresponding region of human hyaluronidase isozyme-1/HYAL1.

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality: Monoclonal

Applications: Research Grade Biosimilar

Target: Hyaluronidase PH-20, Spermatocyte surface protein PH-20, HYAL3, Spermatocyte adhesion molecule 1, PH20, SPAM1, Hyal-PH20, Hyaluronoglucosaminidase PH-20

Purification: Purified by Nickel column.


Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: CAS:2636716-20-4

Note: For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-pan SARS-CoV-2 S2/Spike glycoprotein 2 Antibody (CC99.103)

SDS-PAGE