

Product Information

Catalog No.	BS12139
Product name	Research Grade Gulgafafusp Alfa

Background

Host species:

Isotype: humanglucagonlikepeptide1receptor, GLP1R fused at the N-terminus of each light chain (53-266) via the peptide linker 30GSGSATGGSG SGASSGSGSATGS52 to a glucagon-like peptide-1 (GLP-1) fragment

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality: Monoclonal

Applications: Research Grade Biosimilar

Target: GLP-1 receptor, Glucagon-like peptide 1 receptor, GLP1R, GLP-1-R, GLP-1R, GLP-1(7-36), Incretin hormone, GLP-2, GCG, GLP-1(7-37), OXM, Pro-glucagon, GLP-1, GRPP, OXY

Purification: Purified by Nickel column.


Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: CAS:2642374-02-3

Note: For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Gulgafafusp Alfa.

SDS-PAGE