


Product Information	
Catalog No.	BS12144
Product name	Anti-Human VSIR/VISTA Reference Antibody (BMS767)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	BMS767
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	V-set domain-containing immunoregulatory receptor, C10orf54, VSIR, VISTA, SISP1, Stress-induced secreted protein-1, Platelet receptor Gi24, V-set immunoregulatory receptor, Sisp-1, V-type immunoglobulin domain-containing suppressor of T-cell activation
Expression system	Mammalian Cells
Accession	Q9H7M9
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	V-setdomain-containingimmunoregulatoryreceptor, C10orf54, VSIR, VISTA, SISP1, Stress-inducedsecretedprotein-1, PlateletreceptorGi24, V-setimmunoregulatoryreceptor, Sisp-1, V-typeimmunoglobulindomain-containingsuppressorofT-cellactivation
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human VSIR/VISTA Antibody (BMS767)

SDS-PAGE