

Product Information

Catalog No.	BS12149
Product name	Anti-Human CLDN6 & CD3 Reference Antibody (6PHU3)

Background

Host species:

Isotype: bi-(scFv)₂

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 0.9mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality: Monoclonal

Applications: Research Grade Biosimilar

Target: Claudin-6, CLDN6, Skullin, T3E, T-cellsurfaceantigenT3/Leu-4epsilonchain, CD3e, CD3E, T-cellsurfaceglycoproteinCD3epsilonchain

Purification: Purified by Nickel column.

Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: Bispecific

Note: For research use only. Not suitable for clinical or therapeutic use.