


Product Information	
Catalog No.	BS12184
Product name	Anti-Human SIGLEC8/SAF-2 Reference Antibody (1H10)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1H10
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	SAF2, Sialoadhesin family member 2, SAF-2, Siglec-8, SIGLEC8, Sialic acid-binding Ig-like lectin 8
Expression system	Mammalian Cells
Accession	Q9NYZ4
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human SIGLEC8/SAF-2 Antibody (1H10)

SDS-PAGE