


Product Information	
Catalog No.	BS12207
Product name	Anti-Human CD262/TNFRSF10B/DR5 Reference Antibody (TAS266)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Clonality	Monoclonal
Isotype	VHH-VHH-VHH-VHH
Clone ID	TAS266
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TRAILR2, TRAIL-R2, Death receptor 5, CD262, DR5, Tumor necrosis factor receptor superfamily member 10B, ZTNFR9, TRAIL receptor 2, TRICK2, TNFRSF10B, KILLER, TNF-related apoptosis-inducing ligand receptor 2
Expression system	Mammalian Cells
Accession	O14763
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Givastomig.

SDS-PAGE