

Product Information

Catalog No.	BS12213
Product name	Anti-Human CD20 Reference Antibody (Rituximab, RUO)

Background

Host species: Human

Isotype: IgG1-kappa

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality: Monoclonal

Applications: Research Grade Biosimilar

Target: B-lymphocytesurfaceantigenB1,Membrane-spanning4-domainssubfamilyAmember1,MS4A1,LeukocytesurfaceantigenLeu-16,B-lymphocyteantigenCD20,Bp35,CD20

Purification: Protein A/G purified from cell culture supernatant.


Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names:

Note: For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Lunsekimig.

SDS-PAGE