


Product Information	
Catalog No.	BS12249
Product name	Anti-Human CD161 Reference Antibody (Rezorstobart, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CLEC5B, KLRB1, CD161, HNKRP-1a, C-type lectin domain family 5 member B, NKRP1A, Natural killer cell surface protein P1A, Killer cell lectin-like receptor subfamily B member 1, NKR-P1A
Expression system	Mammalian Cells
Accession	Q12918
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	2762201-86-3
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Rezorstobart

SDS-PAGE