


Product Information	
Catalog No.	BS12256
Product name	Anti-BKPyV VP1 Reference Antibody (Traxivitug, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	BK polyomavirus (BKPyV) (Human polyomavirus 1)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-lambda2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Major capsid protein VP1, Major structural protein VP1
Expression system	Mammalian Cells
Accession	P03088
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	2770852-89-4, MAU868
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Traxivitug.

SDS-PAGE

Traxivitug binds with Major capsid protein VP1

Detects recombinant BKPyV VP1 (Catalog No: VK452011) in indirect ELISA.

