

Product Information

Catalog No.	BS12272
Product name	Research Grade Sifarobocept

Background

Host species:

Isotype: HumanROBO2fragment(1-203), fused to human IgG1Fc fragment

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality:

Applications: Research Grade Biosimilar

Target: Roundabout homolog 2, KIAA1568, ROBO2

Purification: Protein A/G purified from cell culture supernatant.

Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: CAS:2760175-31-1

Note: For research use only. Not suitable for clinical or therapeutic use.