

Product Information	
Catalog No.	BS12297
Product name	Anti-Human IL38/IL1F10 Reference Antibody (IMM-ONC-01)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	IMM-ONC-01
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Interleukin-1 family member 10, IL-1F10, Family of interleukin 1-theta, FIL1 theta, Interleukin-1 HY2, IL-1HY2, Interleukin-1 theta, IL-1 theta, Interleukin-38, IL-38, IL1F10, FIL1T, IL1HY2, IL38, FKSG75, UNQ6119/PRO20041
Expression system	Mammalian Cells
Accession	Q8WWZ1
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Note	For research use only. Not suitable for clinical or therapeutic use.