


Product Information	
Catalog No.	BS12299
Product name	Anti-HBsAg Reference Reference Antibody (OST 577)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Hepatitis B virus genotype A2 subtype adw (isolate Japan/Nishioka/1983) (HBV-A)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-lambda
Clone ID	OST 577
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Large envelope protein, L glycoprotein, L-HBsAg, LHB, Large S protein, Large surface protein, Major surface antigen, S
Expression system	Mammalian Cells
Accession	P03142
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	OST577, OST577, OST-577
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-HBsAg Reference Antibody (OST 577).

SDS-PAGE