

Product Information	
Catalog No.	BS12300
Product name	Anti-Human CD307e/FCRL5 Reference Antibody (DFRF 4539A)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa
Clone ID	DFRF 4539A
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Fc receptor-like protein 5, FcR-like protein 5, FcRL5, BXMAS1, Fc receptor homolog 5, FcRH5, Immune receptor translocation-associated protein 2, CD307e, FCRL5, FCRH5, IRTA2, UNQ503/PRO820
Expression system	Mammalian Cells
Accession	Q96RD9
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	DFRF4539A, DFRF-4539A, DFRF4539A, RG7598, RG-7598, RG7598
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data

Detects recombinant human CD307e/FCRL5 (Catalog No: HV698011) in indirect ELISA

Anti-Human CD307e/FCRL5 Antibody (DFRF 4539A)
binds with CD307e/FCRL5

