

Product Information

Catalog No.	BS12329
Product name	Research Grade Goflikicept

Background

Host species:

Isotype: Interleukin1ReceptorAccessoryProtein(human)fusionproteinwithImmunoglobulinFcFragment.

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality: Monoclonal

Applications: Research Grade Biosimilar

Target: IL1F1,Interleukin-1alpha,IL1A,Hematopoietin-1,IL-1alpha,Interleukin-1beta,IL1B,Catabolin,IL1F2,IL-1beta

Purification: Protein A/G purified from cell culture supernatant.

Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names:

Note: For research use only. Not suitable for clinical or therapeutic use.