


Product Information	
Catalog No.	BS12351
Product name	Anti-Human CCR8 Reference Antibody (Lanerkitug, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-lambda2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CKR-L1, CMKBR8, C-C CKR-8, GPR-CY6, CCR8, CMKBRL2, TER1, CCR-8, Chemokine receptor-like 1, CKRL1, CDw198, GPRCY6, C-C chemokine receptor type 8, CC-CKR-8, CC chemokine receptor CHEMR1
Expression system	Mammalian Cells
Accession	P51685
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	2839417-54-6
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Lanerkitug.

SDS-PAGE