

Product Information

Catalog No.	BS12355
Product name	Research Grade Brivekimig

Background

Host species:

Isotype: ImmunoglobulinGsinglechainof5VH(VH-VH-VH'-VH''-VH')monomer

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality: Monoclonal

Applications: Research Grade Biosimilar

Target: CD252,OX40L,TXGP1,OX40ligand,GlycoproteinGp34,TAXtranscriptionally-activatedglycoprotein1,TNFSF4,Tumornecrosisfactorligandsuperfamilymember4,TNF,Tumornecrosisfactorligandsuperfamilymember2,N-terminalfragment,ICD2,NTF,TNF-a,TNF-alpha,Tumornecrosisfactor,TNFSF2,TNFA,Cachectin,ICD1,ALB,Albumin

Purification: Protein A/G purified from cell culture supernatant.


Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: Trispecific,Pentavalent,CAS:2834733-45-6

Note: For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Brivekimig

SDS-PAGE