


Product Information	
Catalog No.	BS12390
Product name	Anti-Influenza A virus HA/Hemagglutinin Reference Antibody (Sonavibart, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Influenza A virus
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Expression system	Mammalian Cells
Accession	P03455
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	2769699-23-0
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Sonavibart

SDS-PAGE