

Product Information

Catalog No.	BS12408
Product name	Research Grade Palverafusp Alfa

Background

Host species: Humanized

Isotype:

Human VEGFR-1 fragment, fused via a (G4S)₃ peptide linker to the N-terminus of both heavy chains of a humanized IgG1-kappa anti-human PD-L1.

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH 7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality:

Applications: Research Grade Biosimilar

Target: FLT1, Vascular endothelial growth factor receptor 1, VEGFR1, Vascular permeability factor receptor, FRT, Tyrosine-protein kinase receptor FLT, VEGFR-1, FLT, FLT-1, Tyrosine-protein kinase FRT, Fms-like tyrosine kinase 1, B7-H1, Programmed cell death 1 ligand 1, PDCD1 ligand 1, PDCD1L1, B7 homolog 1, PDCD1LG1, PDL1, hPD-L1, Programmed death ligand 1, B7H1, PD-L1, CD274

Purification: Protein A/G purified from cell culture supernatant.


Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: CAS:2850355-87-0

Note: For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Palverafusp Alfa.

SDS-PAGE

Palverafusp Alfa binds with PD-L1

Detects recombinant PDL1 (Catalog No: HV974011) in indirect ELISA.

