


Product Information	
Catalog No.	BS12443
Product name	Anti-Human CD315/PTGFRN Reference Antibody (AG02-ADC)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Isotype	IgG
Clone ID	AG02-ADC
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD9 partner 1, PTGFRN, EWI-F, CD9P-1, CD315, CD9P1, FPRP, Prostaglandin F2 receptor negative regulator, Prostaglandin F2-alpha receptor-associated protein, KIAA1436, Glu-Trp-Ile EWI motif-containing protein F, Prostaglandin F2-alpha receptor regulatory protein, EWIF
Expression system	Mammalian Cells
Accession	Q9P2B2
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	AG02, AG-02, AG02
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human CD315/PTGFRN Antibody (AG02-ADC)

SDS-PAGE