


Product Information	
Catalog No.	BS12488
Product name	Anti-Human MAGEA4 Reference Antibody (Imc-C103C)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	TCR-scFv
Clone ID	Imc-C103C
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Cancer/testis antigen 1.4, MAGEA4, MAGE4, CT1.4, Melanoma-associated antigen 4, MAGE-X2 antigen, MAGE-41 antigen, MAGE-4 antigen
Expression system	Mammalian Cells
Accession	P43358
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	Imc-C103C, ImcC103C, ImcC103C
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-HRSV F/Fusion glycoprotein F0 Antibody (IDB-002)

SDS-PAGE