

Product Information	
Catalog No.	BS12493
Product name	Anti-Human TNFSF15/TL1A Reference Antibody (TEV-574)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Isotype	IgG1-lambda
Clone ID	TEV-574
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TNFSF15, Vascular endothelial cell growth inhibitor, VEGI, TL1, TNF ligand-related molecule 1, Tumor necrosis factor ligand superfamily member 15
Expression system	Mammalian Cells
Accession	O95150
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	TEV-574, TEV574, TEV574, TEV'574
Note	For research use only. Not suitable for clinical or therapeutic use.