

Product Information	
Catalog No.	EL10008
Product name	Lupartumab ELISA Kit
Applications	ELISA
Accession	O95274
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 ug/ml
Range	0.31-5 ug/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	BAY 1112623, BAY 1129980, , 1640971-88-5
Shipping	2-8 deg C
Note	For Research Use Only.