

Product Information	
Catalog No.	EL10030
Product name	Namimumab ELISA Kit
Applications	ELISA
Accession	P04141
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 ug/ml
Range	0.31-5 ug/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	MT203, 1206681-39-1
Shipping	2-8 deg C
Note	For Research Use Only.

Background

Namimumab (alternative identifier MT203) is a human monoclonal antibody (class IgG1 kappa) that targets granulocyte macrophage-colony stimulating factor (GM-CSF), which is also known as colony stimulating factor 2 (CSF2). Namimumab was produced by Micromet Inc and is under development by Takeda Pharmaceuticals International. It is currently being researched for application in rheumatoid arthritis (RA) and psoriatic arthritis. Clinical trials investigating the therapeutic utility of Namimumab have include phase I and phase II clinical trials to establish the safety, tolerability and preliminary therapeutic utility of the antibody in plaque psoriasis and rheumatoid arthritis.