

Product Information	
Catalog No.	EL10035
Product name	Trastuzumab ELISA Kit
Applications	ELISA
Accession	P04626
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	70.12 ng/mL
Range	78.13 - 5, 000 ng/mL
Recovery	80-120%
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20deg C, the rest reagents should be store at 4deg C.
Alternative Names	timigutuzumab, 4D5-8, 4D5V8, Herceptin, rhuMab HER2, 180288-69-1
Shipping	2-8 deg C
Note	For Research Use Only.