

Product Information	
Catalog No.	EL10037
Product name	Zenocutuzumab ELISA Kit
Applications	ELISA
Accession	P04626 & P21860
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 ug/ml
Range	0.31-5 ug/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	MCLA-128, PB4188, R040517, 1969309-56-5
Shipping	2-8 deg C
Note	For Research Use Only.

Background

Zenocutuzumab, an investigational bispecific antibody, has shown early efficacy in an ongoing phase I/II study of patients with tumors harboring NRG1 fusions. The data so far appear particularly promising in NRG1 fusion-positive pancreatic ductal adenocarcinoma, which tends to occur in younger patients.