

Product Information	
Catalog No.	EL10098
Product name	Vofatamab ELISA Kit
Applications	ELISA
Accession	P22607
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 ug/ml
Range	0.31-5 ug/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	MFGR1877S, 1312305-12-6
Shipping	2-8 deg C
Note	For Research Use Only.

Background

Vofatamab is a novel monoclonal antibody specific for fibroblast growth factor receptor 3 (FGFR3) that is being developed to target FGFR3-positive tumors.