

Product Information	
Catalog No.	EL10151
Product name	Tibulizumab ELISA Kit
Applications	ELISA
Accession	Q16552 & Q9Y275
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 ug/ml
Range	0.31-5 ug/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	LY 3090106, LY-3090106, LY3090106, 1849636-24-3
Shipping	2-8 deg C
Note	For Research Use Only.

Background

Tibulizumab (LY-3090106) is a tetravalent bispecific antibody targeting BAFF and IL-17A for the treatment of autoimmune disease.