

Product Information	
Catalog No.	EL10153
Product name	Tiragolumab ELISA Kit
Applications	ELISA
Accession	Q495A1
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	31.25 ng/mL - 2000 ng/mL
Recovery	75-110%
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20deg C, the rest reagents should be store at 4deg C.
Alternative Names	MTIG-7192-A, MTIG-7192A, MTIG7192A, RG-6058, RO-7092284, RO7092284, 1918185-84-8
Shipping	2-8 deg C
Note	For Research Use Only.

Background

Tiragolumab is an investigational, fully human monoclonal IgG1 antibody designed to prevent PVR-TIGIT interaction.