

Product Information	
Catalog No.	EL10162
Product name	Ianalumab ELISA Kit
Applications	ELISA
Accession	Q96RJ3
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 ug/ml
Range	0.31-5 ug/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	VAY-736, 1929549-92-7
Shipping	2-8 deg C
Note	For Research Use Only.

Background

Ianalumab (VAY736) is an anti-B-cell activating factor (BAFF) receptor fully human IgG1 monoclonal antibody, engineered for direct ADCC-mediated B-cell depletion, thus providing a dual mode of action and targeted approach to treat pSS. The current phase 2b study aimed at establishing a dose-response relationship over a wide range of VAY736 doses, using change from baseline (BL) in the EULAR Sjogren's Syndrome Disease Activity index (ESSDAI) over 24 Weeks (Wks) as the primary endpoint. Here we report Wk24 efficacy and safety. The study is ongoing with a second blinded treatment period up to Wk52.