

Product Information	
Catalog No.	EL10167
Product name	Navicixizumab ELISA Kit
Applications	ELISA
Accession	Q9NR61 & P15692
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 ug/ml
Range	0.31-5 ug/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	OMP-305B83, 1638338-43-8
Shipping	2-8 deg C
Note	For Research Use Only.

Background

Navicixizumab is an anti-DLL4/VEGF bispecific antibody that has demonstrated robust in vivo antitumor efficacy in a number of solid tumor xenografts.