

Product Information	
Catalog No.	EL10190
Product name	Monkeypox Virus (MPXV) A35R Antigen ELISA Kit
Applications	ELISA
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	158.4 pg/mL
Range	375 pg/mL - 24, 000 pg/mL
Recovery	70-100%
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20deg C, the rest reagents should be store at 4deg C.
Shipping	2-8 deg C
Note	For Research Use Only.