

Product Information	
Catalog No.	EL10192
Product name	Anti-Eculizumab Neutralizing Antibody ELISA kit
Applications	ELISA
Accession	P01031
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	15.625-1, 000 ng/ml
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	5G1.1, h5G1.1HuG2/G4, 219685-50-4
Shipping	2-8 deg C
Note	For Research Use Only.