

Product Information	
Catalog No.	EL10247
Product name	Anti-Infliximab ELISA Kit
Applications	ELISA
Accession	P01375
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	156.25 - 10, 000 ng/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	TA-650, cA2, 170277-31-3
Shipping	2-8 deg C
Note	For Research Use Only.