

Product Information	
Catalog No.	EL10275
Product name	Anti-HRSV F/Fusion glycoprotein F0 Human IgG ELISA Kit
Applications	ELISA
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.88 ng/mL
Range	1.56 ng/mL - 100 ng/mL
Recovery	80-120%
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20deg C, the rest reagents should be store at 4deg C.
Shipping	2-8 deg C
Note	For Research Use Only.