

Product Information	
Catalog No.	EL10337
Product name	Anti-MARV Envelope glycoprotein (GP) Human IgM ELISA Kit
Applications	ELISA
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	23 ng/mL
Range	93.75 - 6, 000 ng/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Shipping	2-8 deg C
Note	For Research Use Only.