

Product Information	
Catalog No.	EL10377
Product name	SARS-CoV-2 Surrogate Virus Neutralization Test Kit (JN.1)
Applications	ELISA
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	16.81 ng/mL
Range	31.25 - 2000 ng/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Shipping	2-8 deg C
Note	For research use only.