

Product Information	
Catalog No.	EN10000
Product name	CAT/Catalase
Applications	Proteomics
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>90% as determined by SDS-PAGE.
Expression system	E. coli
Protein length	CAT / Catalase is cloned from human cell lines and expressed in E.coli.
Predicted molecular weight	62.07 kDa
Nature	Recombinant
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details.
Species	Human
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.