


| Product Information   |                                                                                                                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | FC10020                                                                                                                                                                                                                                                      |
| Product name          | Anti-Human VSIR/VISTA Antibody (SAA0172)                                                                                                                                                                                                                     |
| Applications          | FCM                                                                                                                                                                                                                                                          |
| Species reactivity    | Human                                                                                                                                                                                                                                                        |
| Host species          | Human                                                                                                                                                                                                                                                        |
| Clonality             | Monoclonal                                                                                                                                                                                                                                                   |
| Isotype               | IgG1, kappa                                                                                                                                                                                                                                                  |
| Clone ID              | SAA0172                                                                                                                                                                                                                                                      |
| Form                  | Liquid                                                                                                                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                                                                                                                           |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                      |
| Purity                | >95% as determined by SDS-PAGE.                                                                                                                                                                                                                              |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                                                                                                                          |
| Target                | V-set domain-containing immunoregulatory receptor, C10orf54, VSIR, VISTA, SISP1, Stress-induced secreted protein-1, Platelet receptor Gi24, V-set immunoregulatory receptor, Sisp-1, V-type immunoglobulin domain-containing suppressor of T-cell activation |
| Accession             | Q9H7M9                                                                                                                                                                                                                                                       |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                                                                                                                            |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.                                                                                                 |
| Note                  | For research use only.                                                                                                                                                                                                                                       |

### Image Data



SDS-PAGE for Anti-Human VSIR/VISTA Antibody (SAA0172)

SDS-PAGE