

Product Information	
Catalog No.	FC10041
Product name	Anti-Human SLC1A5 Antibody (SAA0181)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0181
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.58 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Baboon M7 virus receptor, SLC1A5, Solute carrier family 1 member 5, ASCT2, RDR, RD114/simian type D retrovirus receptor, Sodium-dependent neutral amino acid transporter type 2, ATB(0), Neutral amino acid transporter B(0), M7V1, RDRC
Accession	Q15758
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.