


Product Information	
Catalog No.	FC10043
Product name	Anti-Human CD119/IFNGR1 Antibody (SAA0044)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0044
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IFN-gamma-R1, Interferon gamma receptor 1, CDw119, Interferon gamma receptor alpha-chain, IFN-gamma receptor 1, IFNGR1, IFN-gamma-R-alpha, CD119
Accession	P15260
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD119/IFNGR1 Antibody (SAA0044)

SDS-PAGE