


Product Information	
Catalog No.	FC10055
Product name	Anti-Human PCSK9 Antibody (SAA0158)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG2, lambda
Clone ID	SAA0158
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.18 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Proprotein convertase subtilisin/kexin type 9, Subtilisin/kexin-like protease PC9, NARC-1, PC9, Proprotein convertase 9, Neural apoptosis-regulated convertase 1, NARC1, PCSK9
Accession	Q8NBP7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human PCSK9 Antibody (SAA0158)

SDS-PAGE