


Product Information	
Catalog No.	FC10062
Product name	Anti-Human MUC16/CA125 Antibody (3D5)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	3D5
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Mucin-16, CA-125, Ovarian carcinoma antigen CA125, CA125, MUC-16, Ovarian cancer-related tumor marker CA125, MUC16
Accession	Q8WXI7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human MUC16/CA125 Antibody (3D5)

SDS-PAGE