


Product Information	
Catalog No.	FC10082
Product name	Anti-Human IL1RL1/ST2 Antibody (SAA0139)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0139
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	T1, DER4, Interleukin-1 receptor-like 1, ST2, IL1RL1, Protein ST2
Accession	Q01638
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human IL1RL1/ST2 Antibody (SAA0139)

SDS-PAGE