

Product Information	
Catalog No.	FC10094
Product name	Anti-Human HLA-A, B, C Antibody (W6/32)
Applications	ELISA, FCM, ICC, IHC-Fr, IP, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	W6/32
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	6 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.