


Product Information	
Catalog No.	FC10099
Product name	Anti-Human GPC3 Antibody (SAA0134)
Applications	ELISA, FCM, WB
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0134
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	GPC3, MXR7, GTR2-2, Glypican-3, Intestinal protein OCI-5, OCI5
Accession	P51654
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human GPC3 Antibody (SAA0134)

SDS-PAGE