

| Product Information   |                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | FC10116                                                                                                                                                      |
| Product name          | Anti-Human ERBB3/HER3 Antibody (SAA0122)                                                                                                                     |
| Applications          | FCM                                                                                                                                                          |
| Species reactivity    | Human                                                                                                                                                        |
| Host species          | Human                                                                                                                                                        |
| Clonality             | Monoclonal                                                                                                                                                   |
| Isotype               | IgG1, kappa                                                                                                                                                  |
| Clone ID              | SAA0122                                                                                                                                                      |
| Form                  | Liquid                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                           |
| Concentration         | 1 mg/ml                                                                                                                                                      |
| Purity                | >95% as determined by SDS-PAGE.                                                                                                                              |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                          |
| Target                | ERBB3, HER3, Proto-oncogene-like protein c-ErbB-3, Tyrosine kinase-type cell surface receptor<br>HER3, Receptor tyrosine-protein kinase erbB-3               |
| Accession             | P21860                                                                                                                                                       |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                            |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| Note                  | For research use only.                                                                                                                                       |